


LQT6280LO 0459

IN-LINE RED CELL FILTRATION



BLOOD PROCESSING



TECHNICAL SPECIFICATIONS

| Bag(s) | Filter(s) | Collection volume | Packaging | |
|------------|----------------|-------------------|----------------------------|------------------|
| 4 PVC-DEHP | LEUCOFLEX LCRD | 450 mL +/- 10 % | 1 device / peelable sachet | 12 devices / box |

Included items

- Needle (16G)
- Bactivam, Secuvam and Vacuam
- Break-Away Cannula in the Sampling Line

Solutions

- 63 mL CPD
- 100 mL SAGM

Sterilisation

- Steam

Labels

- English
- ISBT128

Shelf life

- 2 years

Instructions for use

NOTEURO04

Final products

- Plasma
- Buffy-Coat
- Leucoreduced Red Cell Concentrate

Performances

- AABB guideline
- European guideline
- Processing up to 24 hours at room temperature